

### Care kit for Valplast

*Everything you need for easy cleaning or your valplast flexible denture.....£58*

The care kit has been put together to make it easy and more convenient to look after a valplast flexible denture. It includes the items required to keep the denture looking and feeling like new. Supplied in a convenient zip top, re-usable presentation kit bag.

Content:

- 1x Val-clean - (approx 12 months supply)
- 1x Sonic Denture Cleaner
- 1x Cleaning solution store
- 1x Second generation Flexible Denture brush with Travel Tube
- 1x Presentation Care Kit Bag



### Val-Clean Concentrated Denture Cleaner

*Box of 12 monthly sachets, concentrated power - approx 12 months supply .....£29*

Val-Clean is a non-abrasive cleaner specifically designed for the hygienic cleaning and maintenance of valplast appliances. Val-Clean is supplied in boxes of 12 sachets, each sachet will make four weeks supply of solution. A mixing spoon is included with each box, one spoonful of Val-Clean is mixed with approximately 7oz (250ml) of warm water - this solution will last for 7- 10 days before needing to be replaced. Val-Clean is designed specifically for cleaning Valplast dentures. A daily 15 minute soaking has proven to be effective in maintaining the appearance and cleanliness of Valplast. Overnight soaking of Val-Clean solution will not harm the appliance.

### Sonic Denture Cleaner

*The way to clean your Valplast flexible partial.....£29*

Through the chemical action of Val-Clean or FDC Flexible Denture Cleaner, plus sonic vibration from the Sonic Denture Cleaner, tartar, stain, film, bacteria, and food particles are quickly and easily removed from your Valplast denture. Simply rinse your denture under running water before immersing it in the cleaning solution in your Sonic Cleaner. Switch on the sonic Cleaner for 10- 15 minutes. After cleaning the denture can remain soaking in the cleaning solution until required for use.

Turn over for more products...

## Flexible Denture Brush

*Specially developed to be used on flexible partial dentures.....£6*

Flexible Denture Brush - Specially developed to be used on flexible partial dentures. The bristles are made from 100% highest quality silicone, soft enough not to damage your flexible partial denture but firm enough to remove food particles. It is supplied in a travel tube for the safe transport and storage of the brush. (Note: Traditional toothbrushes may remove the polish and roughen the surface of flexible partial dentures in time)

## FDC Flexible Denture Cleaner and Disinfectant

*Box of 12 Weekly Sachets, concentrated power - 3 months supply (approx).....£12*

FDC Flexible Denture Cleaner and Disinfectant for use with Flexible Partial Dentures, Acrylic Dentures, Partial Dentures, Night Guards and Retainers. Not only does FDC clean but it also disinfects as it cleans, mint flavoured for an added fresh feel to your denture. While disinfecting, FDC will also remove coffee, tea, tobacco, fruit and iron stains, calcium, calculus, tartar and plaque. The result - clean and disinfected dentures. According to several studies conducted in universities and hospitals 65% of denture wearers suffer from problems caused by Candida Albicans and other pathogenic micro-organisms present on removable appliances. These germs can cause serious oral inflammation. Unless eliminated, this condition can lead to denture intolerance. Two independent laboratory studies confirm that FDC Flexible Denture Cleaner totally eliminates Candida Albicans, Actinomyces Viscosus, Streptococcus mutans and streptococcus pyogenes from the denture in just 20 minutes.

## Flexible Denture Cleaning Container

*Robust cleaning container, with a removable inner basket.....£5*

Designed for the safe storage and hygienic cleaning of Valplast Flexible Dentures. Mix either Val-Clean Concentrated Denture Cleaner or FDC Flexible denture cleaner according to the instructions on the packet. Pour approximately 200ml of the prepared solution into the cleaning container and totally immerse the denture for about 10 to 15 minutes to achieve a clean fresh denture. The cleaning container is a safe and secure storage place for your denture when overnight soaking is used. The removable basket can be used to remove the denture from the solution and to rinse it before replacing it in your mouth. The Cleaning Container should be thoroughly cleaned before replacing the cleaning solution.

## Cleaner Solution Store

*Safe storage of the prepared Val-Clean & FDC (100ml).....£2*

The store is ideal for use when traveling and to avoid accidental spillage during the cleaning and disinfecting process. Val-Clean & FDC remains active for 7-10 days making it possible to use the solution daily, the store makes an ideal place to keep the solution. The Cleaner Solution Store has been designed to take most sizes of Valplast Partial Denture complete with 100ml of prepared solution.

## Care Kit Bag for Valplast

*Zip top, re-usable care kit bag to conveniently keep all of your Valplast items together when not in use .....£4*