




Variolink<sup>®</sup> Veneer  
ESTHETIC RESIN CEMENT

ivoclar  
vivadent<sup>®</sup>

# EXCEPTIONAL ESTHETICS

## Value Shading

Variolink Veneer offers seven light cure value shades to enhance the esthetics of the final restoration. The Medium Value (MV 0) is very translucent and has little effect on the final restoration. The High Values allow the clinician to brighten or increase the value and the Low Values provide a "warming" effect to the final restoration.



## Variolink Veneer Light Cure

All Variolink Veneer are amine-reduced providing a high degree of color stability ensuring outstanding esthetics.



Variolink<sup>®</sup> Veneer

# EXTRAORDINARY CLEAN UP

## Clean Up

Variolink Veneer offers the clinician a fast, easy solution to clean up excess cement during the adhesive bonding procedure. No more grinding, burs, or tedious clean up required. Simply load the restorations, seat, cure, and peel away the excess. Variolink Veneer provides the unique combination of the security of adhesive bonding and the easy clean up of a conventional cement.




Clinical photography courtesy of Dr. David Hornbrook

## Physical Properties


Variolink Veneer's microfill chemistry provides key benefits including low wear rates, high translucency, high polishability and strong bonds. Variolink Veneer demonstrates significantly lower wear values compared to competitive resin cements currently on the market ensuring the long term clinical success of your esthetic restorations.

### TRANSLUCENCY



\* No registered trademarks of Ivoclar Vivadent, Inc. 2003, R&D Ivoclar Vivadent, Inc.

### WEAR



\* No registered trademarks of Ivoclar Vivadent, Inc.

UAB School of Dentistry, October 2003  
Shiro Suzuki, DDS, PhD, FAD  
Professor of Dentistry  
Director of Clinical Research  
Department of Prosthodontics and Biomaterials

# Variolink® Veneer

## DELIVERY FORMS

### Assortment/Syntac



- 7 syringes Variolink Veneer, 1.5 ml each, in 7 shades
- 7 syringes Variolink Veneer Try-In, 1.5 ml each, in 7 shades
- 1 bottle Monobond S, 5 g
- 1 syringe Liquid Strip, 2.5 g
- 1 syringe Total Etch, 2 g
- 1 bottle Syntac Primer, 3 g
- 1 bottle Syntac Adhesive, 3 g
- 1 bottle Heliobond, 6 g
- Various accessories

### Assortment/Excite



- 7 syringes Variolink Veneer, 1.5 ml each, in 7 shades
- 7 syringes Variolink Veneer Try-In, 1.5 ml each, in 7 shades
- 1 bottle Monobond S, 5 g
- 1 syringe Liquid Strip, 2.5 g
- 1 syringe Total Etch, 2 g
- 12 Excite single dose vessels, 0.1 g each
- 1 syringe Total Etch, 2 g
- Various accessories

### Refills/Variolink Veneer



- 1 syringe Variolink Veneer, 1.5 ml, available in the following shades:

High Value	+3
High Value	+2
High Value	+1
Medium Value	0
Low Value	-1
Low Value	-2
Low Value	-3

### Refills/Variolink Veneer Try-In



- 1 syringe Variolink Veneer Try-In, 1.5 ml, available in the following shades:

High Value	+3
High Value	+2
High Value	+1
Medium Value	0
Low Value	-1
Low Value	-2
Low Value	-3

Descriptions and data constitute no warranty of attributes.

Printed in Liechtenstein  
© Ivoclar Vivadent AG, Schaan  
Liechtenstein  
598376/0905/e/BVD

Ivoclar Vivadent AG Clinical  
Bendererstr. 2  
FL-9494 Schaan  
Fürstentum Liechtenstein  
Tel. +423 / 235 35 35  
Fax +423 / 235 33 60  
www.ivoclarvivadent.com

ivoclar  
vivadent®